



# **MANTECH SPOTLIGHT**

A complete system for only

**\$22,995 (Optical Probe)**

**\$23,995 (Membrane Probe)**

Valid until June 30, 2016

## **MANTECH PC-BOD Mini Read Only System**

Capable of Automating the BOD5 and BOD7 methods.

System includes Automax73 Sampler, PC controller, Software, YSI Optical or Membrane Dissolved Oxygen probe, installation and user training.

**REQUEST A QUOTE**  
**1-866-763-2122**



- Specifically designed for laboratories analyzing a small number of samples per day
- Use glass or plastic 300 mL or glass 60 mL BOD bottles
- Analyze BOD and CBOD in the same run
- Easy to use PC-BOD software
- Can automatically add seed, inhibitor and dilution water as required with optional pumps

